

EPI_REST

PATIENT POSITION	
Patient Entry	Head First
Patient Position	Supine
Coil Configuration	8HRBRAIN
Plane	AXIAL
Series Description	EPI_REST
SCAN TIMING	
Flip Angle	90
TE	30.0
Number of Echoes	1
TR	2000.0
Number of Shots	1
IMAGE ENHANCE	
Filter Choice	None
GATING/TRIGGER	
Auto Trigger Type	Off
FMRI	
Brain Wave Real Time	0
Paradigm String	Research Mode
Paradigm UID	1.2.840.113819.3.20060302.1
Initial State	Control
PSD Trigger	Internal
Slice Order	Interleaved
View Order	Bottom/Up
# of Repetitions REST	1
# of Repetitions ACTIVE	1
# of Dummy Acquisition	5
SAT	
Tag Type	None
TRICKS	
Pause On/Off	On
Auto Subtract	0
Auto SCIC	Off

EPI_REST

IMAGING PARAMETERS	
Imaging Mode	2D
Pulse Sequence	Gradient Echo
Imaging Options	EPI, Asset, FMRI
Acceleration Factor	1.00
SCAN RANGE	
FOV	22.0
Slice Thickness	3.4
Slice Spacing	0.0
ACQ TIMING	
Freq	64
Phase	64
Freq DIR	R/L
NEX	1.00
Phase FOV	1.00
Auto Shim	Auto
Phase Correction	Yes
USER CVS	
User CV0	1.00
MULTI-PHASE	
Slice per Location	180
Delay after Acquisition	0
# of Acquisition	0
Seperate Series	0
Delay after Acquisition without AV	0
Mask Phase	0
Mask Pause	0
ASSET	
Slice Acceleration Factor	1.00
Phase Acceleration Factor	2.00
CONTRAST	
Contrast Yes/No	No

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	Head First	Imaging Mode	2D
Patient Position	Supine	Pulse Sequence	Gradient Echo
Coil Configuration	8HRBRAIN	Imaging Options	None
Plane	AXIAL	PSD Name	grass_ucsd
Series Description	fm_grass_64x64	Acceleration Factor	1.00
SCAN TIMING		SCAN RANGE	
Flip Angle	45	FOV	22.0
TE	6.5	Slice Thickness	3.4
Number of Echoes	1	Slice Spacing	0.0
TR	614.0		
Receiver Bandwidth	31.25		
IMAGE ENHANCE		ACQ TIMING	
Filter Choice	None	Freq	64
		Phase	64
		Freq DIR	R/L
		Auto Shim	Off
		Phase Correction	No
GATING/TRIGGER		USER CVS	
Auto Trigger Type	Off	User CV0	1.00
		User CV1	2.00
FMRI		MULTI-PHASE	
PSD Trigger	Internal	# of Acquisition	0
Slice Order	Interleaved	Seperate Series	0
View Order	Bottom/Up	Mask Phase	0
# of Repetitions REST	0	Mask Pause	0
# of Repetitions ACTIVE	0		
SAT		ASSET	
Tag Type	None	Slice Acceleration Factor	1.00
		Phase Acceleration Factor	1.00
TRICKS		CONTRAST	
Pause On/Off	On	Contrast Yes/No	No
Auto Subtract	0		
Auto SCIC	Off		